

## mikrocid® wipes

Safety on all surfaces.



### mikrocid® AF wipes

The classic choice among rapid disinfection preparations. The alcohol-based wipes are effective within 1 minute against most relevant microbes. The alcohol film dries very quickly and leaves no residues behind. Material compatibility is good, however stress cracks may form in plastics.

**Shelf life\*:** 3 months | **Active ingredients:** Ethanol, Propan-1-ol

Article	Article number	Ø Wipe size	Wipe pieces	Grammage	Surface coverage Ø	Efficacy
Tub	109203	Ø 14 x 18 cm	150	24 g/m <sup>2</sup>	0.3 m <sup>2</sup>	<div>limited</div> <div>BAC VIR TB</div> <div>ADENO NORO FUNG</div> <div>VIRU ZID</div>
Tub Refill	70002062	Ø 14 x 18 cm	150	24 g/m <sup>2</sup>	0.3 m <sup>2</sup>	
Jumbo Tub	70003093	Ø 20 x 20 cm	220	24 g/m <sup>2</sup>	0.5 m <sup>2</sup>	
Jumbo Tub Refill	70003094	Ø 20 x 20 cm	220	24 g/m <sup>2</sup>	0.5 m <sup>2</sup>	
Softpack	70000808	Ø 20 x 20 cm	50	50 g/m <sup>2</sup>	0.5 m <sup>2</sup>	

### mikrocid® sensitive wipes

The specialist for sensitive surfaces. There are many sensitive materials in a hospital and for this reason, mikrocid® sensitive wipes combine a high degree of efficacy with excellent material compatibility. Even ultrasound heads\*\* can be disinfected with mikrocid® sensitive wipes.

**Shelf life\*:** 3 months | **Active ingredient:** QAC  
perfumefree

Article	Article number	Wipe size	Wipe pieces	Grammage	Surface coverage Ø	Efficacy
Tub	70003103	Ø 14 x 18 cm	120	20 g/m <sup>2</sup>	0.3 m <sup>2</sup>	<div>limited</div> <div>BAC VIR LEV</div> <div>NORO</div>
Jumbo Tub	70003104	Ø 20 x 20 cm	200	20 g/m <sup>2</sup>	0.5 m <sup>2</sup>	
Jumbo Tub Refill	70003109	Ø 20 x 20 cm	200	20 g/m <sup>2</sup>	0.5 m <sup>2</sup>	
Softpack	70003101	Ø 20 x 20 cm	50	50 g/m <sup>2</sup>	1 m <sup>2</sup>	
Softpack	70003108	Ø 20 x 20 cm	100	50 g/m <sup>2</sup>	1 m <sup>2</sup>	
Maxi-Softpack	70003102	Ø 25 x 25 cm	80	50g/m <sup>2</sup>	1.5m <sup>2</sup>	

### mikrocid® universal wipes

The all-round wipes – also for disinfecting sensitive communication devices. Smartphones and tablets have nearly constant contact with hands. They can be safely treated with mikrocid® universal wipes without damage to the materials.

**Shelf life\*:** 1 month | **Active ingredients:** Ethanol, Propan-2-ol  
perfumefree

Article	Article number	Wipe size	Wipe pieces	Grammage	Surface coverage Ø	Efficacy
Softpack	70000038	Ø 20 x 20 cm	100	50 g/m <sup>2</sup>	1 m <sup>2</sup>	<div>limited</div> <div>BAC VIR TB</div> <div>LEV ADENO NORO</div>
Maxi-Softpack	70001061	Ø 25 x 25 cm	80	50g/m <sup>2</sup>	1.5m <sup>2</sup>	

### mikrocid® PAA wipes

The powerhouse for cases of outbreak. With the unique oxidising power of peracetic acid, mikrocid® PAA wipes have a sporicidal effect and are also effective against C. difficile in just 5 minutes.

**Shelf life\*:** 28 days | **Active ingredient:** Peracetic acid  
perfumefree

Article	Article number	Wipe size	Wipe pieces	Grammage	Surface coverage Ø	Efficacy
Tub	70003205	Ø 20 x 30 cm	50	40 g/m <sup>2</sup>	3.5 m <sup>2</sup>	<div>BAC VIR TB</div> <div>ADENO NORO SPOR</div> <div>FUN</div>

Use disinfectants safely. Always read the label and product information before use.

\* The shelf life after the container is opened is marked on articles with the following symbol:



\*\* Not suitable for final disinfection of semicritical medical devices.







## schülke wipes

fast.safe.practical.easy

This is how to use the schülke wipes dispenser in the right way:



### Wet wipe dispenser

		schülke wipes mini dispenser bucket (40 wipes dispenser)	schülke wipes 70 dispenser box (70 wipes dispenser)	schülke wipes dispenser bucket (100 wipes dispenser)	schülke wipes 100/111 dispenser box (100 / 111wipes dispenser)	schülke wipes safe & easy bag-in-box (111 wipes in the bag)	schülke wipes safe & easy bagless (130 wipes inclusive)	
		No. 119530 (1 pcs.) No. 70000163 (5 pcs.)	No. 119540 (1 pcs.) No. 70000160 (5 pcs.)	No.119523 (1 pcs.) No. 70000160 (5 pcs.)	No. 119541 (1 pcs.) No. 70000161 (5 pcs.)	No. 119542 (10 pcs.)	No. 70003177 (6 pcs.)	
Picture								
Area		low risk	low risk	low risk	low risk	high risk	high risk	
Operation		Wet wipe dispenser systems for the use in combination with the schülke surface disinfectants for all areas in hospitals, old peoples home, rescue service and doctors practice.					Due to the very easy handling of the single-use wet wipe dispenser systems and the protective bag or bucket the schülke wipes safe & easy are the ideal disinfectant wipes for use in areas with very high risk of contamination with pathogen germs.	
Application		Prepare the system after use or ending of the endurance					Dispose bag and lid after use or ending of the endurance	Dispose bucket after use or ending of the endurance
Amount of wipes / application	schülke wipes mini <b>40 wipes</b> No. 119532	<b>40 wipes</b> 30 x 30 cm						
	schülke wipes <b>70 wipes</b> No. 119505	<b>70 wipes</b> 28 x 23 cm	<b>70 wipes</b> 28 x 23 cm					
	schülke wipes <b>100 wipes</b> No. 119503			<b>100 wipes</b> 30 x 30 cm	<b>100 wipes</b> 30 x 30 cm			
	schülke wipes <b>111 wipes</b> No. 119506			<b>111 wipes</b> 30 x 24 cm	<b>111 wipes</b> 30 x 24 cm			
	schülke wipes <b>safe &amp; easy bag</b> No. 119544					<b>111 wipes (in the bag)</b> 30 x 24 cm		
	<b>130 wipes</b> inclusive in the system						<b>130 wipes (inclusive)</b> 30 x 19 cm	
Solution ready-to-use product / diluted concentrate solution		1,5 litres		2 litres ready-to-use product / 3 litres diluted concentrate solution			2 litres ready-to-use product / 2.5 litres diluted concentrate solution	
Suited for ready-to-use products*:		• antifeck® N liquid • buraton® rapid • mikrozid® AF liquid • mikrozid® universal liquid • mikrozid® sensitive liquid • perform® advanced alcohol EP • perform® classic alcohol EP • terralin® liquid						
Diluted concentrate solution*:		• antifeck® AF (N) • antifeck® extra • perform® classic concentrate Q-Plus • pursept® AF • terralin® protect • terralin® PAA** • TPH protect (Information on other products available on request from our Sales Service department (Phone.: +49 40 521 00-666   E-mail: info@schuelke.com) or from your area sales manager.)						

\* Maximal endurance of 28 days.

\*\* Maximal endurance of up to 12 hours in the schülke wipes mini dispenser bucket.